## ScanAfia XS and XTK Hospital Bed

Safe choice for hospital wards

The quality ScanAfia XS and XTK hospital bed is the safe choice for hospital wards. The bed has been designed for use by healthcare professionals in hospitals, healthcare centres and similar facilities, and it fulfils the requirements of the new EN 60601-2-52:2010 standard for medical beds.

Although it can take loads up to 260 kg, the ScanAfia XS and XTK hospital bed is still light and easy to manoeuvre, thanks to its double castors.

The design places emphasis on both ergonomic working and patient comfort and accessibility; its exceptionally low minimum height makes work easier for nurses and makes getting out of bed easier for patients. The back section slide mechanism also improves patient comfort. The XS bed also features an extremely broad Trendelenburg / anti-Trendelenburg adjustment range.

"Model with comprehensive easy-to-use features."



canAfia X-series	ScanAfia XTK-280/290	ScanAfia XTK-480/490	ScanAfia XS-280/290	ScanAfia XS-480/490
ength, cm	215	215	215	215
Vidth (incl. side rails)	89,5/99,5 cm	89,5/99,5 cm	89,5/99,5 cm	89,5/99,5 cm
Height adjustment (*depending on the castor, pls see matrix below)	3082	3082	3082	3082
Veight, kg	102/104	107/109	104/106	109/111
Recommended mattress size, cm	205 x 78/88 x 13	205 x 78/88 x 13	205 x 78/88 x 13	205 x 78/88 x 13
Safe working load (SWL), kg	260 kg	260 kg	260 kg	260 kg
Central locking standard Ø125mm castors	X	X	Х	Х
Central locking design twin Ø125mm castors	•	•	•	•
Centrally locking standard Ø150mm castors	•	•	•	•
Centrally locking design twin Ø150mm castors	•	•	•	•
Centrally locking twin Ø125mm castors, electric	-	-	•	•
Electric system, 24 V	X	X	X	X
Back section adjustment, 0° - 70°	X	X	X	X
high section adjustment, 0° - 30° (only 4-section models)	-	X	-	X
Calf section adjustment, 0° - 11° (only for 4-part models)	-	X	-	X
rendelenburg/ Anti-Trendelengburg, 12°/20° with gas spring	-	-	X	X
Bed lenghtening of 15cm	•	•	•	•
ifting pole	•	•	•	•
ift support	•	•	•	•
Jrine bag holder	•	•	•	•
ith castor	•	•	•	•
Bedding stand	•	•	•	•
Nonitor stand	•	•	•	•
Nurse locking panel	•	•	•	•
Nurse control panel	•	•	•	•
CU control panel	-	•	-	•
Neal tray	•	•	•	•
Handset	X	X	X	X
Manual emergency lowering of back section	X	X	X	X
Manual (quick) CPR function	-	-	•	•
Back section slide	X	X	X	X
Breathable ABS plastic base with mattress holders	X	X	X	X
Bumper rollers	X	X	X	X
Battery use	-	-	•	•
Side rails, chromed, normal or long version	•	•	•	•
Side rails, painted, normal or long version	•	•	•	•
Bed ends, in different colors	•	•	•	•
DE-marked	X	X	X	X
		.,		
Warranty 2 years	X	X	Χ	X

*I loight adius	atmost rongo	Min 20 am may 20 am (depending on the context are and time)
neight adjus	stment range	Min. 32 cm - max. 82 cm (depending on the castor size and type)

Standard feature = X, Optional feature = X, Not available = X

