

Capre Examination Tables

The Lojer Capre range of examination tables is in a class of its own when it comes to design and function. A large number of examination tables are placed in examination rooms, with no need for manoeuvring – large wheels are therefore unnecessary.

The Lojer Capre E1 examination table is designed for general use and features a one-piece table top. Its exceptionally low minimum height and quiet motor add to its appeal in small examination room environments, such as workplace healthcare.

The Lojer Capre E2 examination table is designed for general use and features a two-piece table top.

The Lojer Capre EG examination table is designed for general gynaecological use and features a two-piece table top.

The product features heel supports and a bowl as standard.

Note! Pictures show optional 360° height adjustment foot bar and castor accessories.

"Best in hi-low range – probably one of the sturdiest range of examination tables in the market".

Capre EG Gynaecological Examination Table with knee supports, castors and 360° height adjustment foot bar



Capre E2 Examination Table with castors and 360° height adjustment foot bar



Capre E1 Examination Table with castors and 360° height adjustment foot bar



Capre Examination Table	Capre E1	Capre E2	Capre EG
Width	67 cm	67 cm	67 cm
Length	196 cm	196 cm	196 cm
Weight	95 kg	96 kg	99 kg
Height adjustment range	45 - 92 cm	45 - 92 cm	45 - 92 cm
Back section length	-	68 cm	68 cm
Back section adjustment, electric	-	0..78°	0..78°
Safe working load (SWL)	210 kg	210 kg	210 kg
Electric system, 24V	X	X	X
Adjustable heel supports	-	-	X
Stainless steel bowl	-	-	X
Foot pedal for height adjustment	-	•	•
Hand set for height adjustment	•	-	-
Special width 80 cm	•	•	-
Battery use	•	•	•
360° height adjustment foot bar	•	•	•
Paper roll holder	•	•	•
Paper cutter	•	•	•
Neck pillow	•	•	•
Side rail	•	•	•
Arm support	•	•	•
Accessory rail	•	•	•
Place for IV stand	•	•	•
IV pole	•	•	•
SKAI Fashion upholstery (13 options)	X	X	X
SKAI Clinica upholstery (3 options)	•	•	•
SKAI Conductive upholstery (1 option)	•	•	•
Castors ø75mm	•	•	•
CE-marking	X	X	X
Warranty 2 years	X	X	X
Country of origin	Finland	Finland	Finland

Standard feature = X, Optional feature = •, Not available = -


Accessories



360° height adjustment foot bar



IV stand



xxxxxxxxxxxxx