


Salli Hydraulic Hospital Beds


Sturdy, functional and mobile

Hydraulic hospital bed Salli is a perfect solution for clinics who ask for sturdy, functional and mobile hospital beds yet of reasonable price. Based on design of popular ScanAfia range Salli hydraulic beds have modular construction thus manage to create models to meet any demands of the medical staff.

*"Functional and
Easy to Move*



Salli H280 Hydraulic Hospital Bed



Salli H480 Hydraulic Hospital Bed

LOJER® For easy care

■ Healthcare | Senior Care | Medical Wellness

SALLI hydraulic beds	H280	H480	H290	H490
Number of sections	2	3	2	3
Width (incl. siderails)	90 cm	90 cm	100 cm	100 cm
Recommended mattress size, cm	205 x 78	205 x 78	205 x 88	205 x 88
Full length	222 cm	222 cm	222 cm	222 cm
Weight	111 kg	114 kg	112 kg	115 kg
Hydraulic height adjustment range	46... 81 cm	46... 81 cm	46... 81 cm	46... 81 cm
Safe Working Load (SWL)	170 kg	170 kg	170 kg	170 kg
Mechanical back section adjustment 0° -70°	X	X	X	X
Gas spring assisted adjustment for back section 0° -70°	•	•	•	•
Mechanical adjustment for knee break 0° -30°	-	X	-	X
Gas spring assisted adjustment for knee break 0° -30°	-	•	-	•
Mechanical foot section adjustment 0° -15°	-	X	-	X
Mechanical lower body part adjustment 0°-15°	•	=	•	=
Centrally locking castors Ø125	X	X	X	X
Centrally locking castors Ø150	•	•	•	•
Individually locking castors Ø125, price reduction	•	•	•	•
One metal bed end, curved/ straight, epoxy-coated/ chromium plated	X	X	X	X
Additional metal bed end, curved/ straight, epoxy-coated/ chromium plated	•	•	•	•
Metal bed end, with tilting mechanism	•	•	•	•
Trendelenburg/Anti-trendelenburg 12°/14°	•	•	•	•
IV-pole straight or curved model	•	•	•	•
Lifting pole	•	•	•	•
Lift support	•	•	•	•
Fifth castor for improved mobility	•	•	•	•
Pair of metal side rails, epoxy-coated/chromium plated	•	•	•	•
Pair of metal side rails, economy model (drop down model)	•	•	•	•
Bumpers	X	X	X	X
CE-marking	X	X	X	X
Warranty	2/10 years	2/10 years	2/10 years	2/10 years
Country of origin	Finland	Finland	Finland	Finland

X = Standard - = not available • = Optional


Accessories



Lifting pole



Side rails



IV-pole