

Transport stretchers

A Reliable Solution for Safe Patient Transport

Lojer transport stretchers are an easy to use, robust and reliable solution for moving patients safely and efficiently. They can be fitted with a number of accessories, for example for first aid use.

We give our products the best warranty on the market, 2/10 years.

"Easy and agile to steer, light construction, excellent price/quality ratio."



Transport Stretchers Lojer 4315


LOJER® For easy care

■ Healthcare | Senior Care | Medical Wellness

Transportation stretcher	Model 4315
Height adjustment, hydraulic	55 / 85 cm
Width	lying surface 63 cm, total width 74 cm
Width, lying surface 70 cm (total width 79 cm)	•
Length	202 cm
Weight, standard model / 70 cm wide model	78 kg / 81 kg
Back section length	72 cm
Back section adjustment with gas spring	X
Trendelenburg, 16°	X
SKAI Fashion upholstery (13 options)	X
Centrally lockable castors Ø125 mm	X
Centrally lockable castors Ø150 mm	•
X-ray pocket under the back section	•
Paper roll holder	•
Oxygen bottle holder for 5 kg bottle	•
Accessory basket for the lower frame	•
Attachment belts	•
Pair of side rails	•
IV stand, straight model	•
Transport stretcher padding 7 cm	•
Transport stretcher padding, anti-static 7 cm	•
Safe Working Load (SWL), 170 kg	X
CE-marking	X
Warranty 2/10 years	X
Country of Origin	Finland

Standard feature = X, Optional feature = •, Not available = -

Accessories



Pair of side rails



IV stand, straight model



Accessory basket for the lower frame