



SCANDIATM 310H

OPERATING TABLE



Scandia is a modern hydraulic operating table with highest functionality and versatility for universal surgical procedures. Scandia 310H provides you functions that have not seen in hydraulic tables before.

- Ergonomic Scandinavian design with lowest height only 670 mm
- Foot controlled hydraulic functions

- with integrated movement selector in the handle
- Gas spring assisted back and leg sections
- Four hydraulic floor locks placed in outer corners for stable positioning during operations
- Easily movable with high patient weight capacity
- X-ray cassette tunnel with optional longitudinal sliding table top for maximum C-arm access
- Wide range of accessories



USING A HYDRAULIC TABLE IS NOW EASIER THAN EVER BEFORE.



SCANDIATM 310H

OPERATING TABLE

SCANDIA 310H FOR GENERAL SURGERY

Four section table top incl. narrow head section, back section and leg section (1-piece).

MATERIALS

Table top and pedals stainless steel and high pressure laminate, base and column casing acryl coated ABS.

OPTIONS

- Sliding top 250 mm, mechanical
- Split leg sections
- Stainless steel covers
- Moulded and viscoelastic mattresses

DIMENSIONS / WEIGHTS

- Overall length with head plate 2170 mm
- Width 550 mm (615 mm over the siderails)

- Patient weight capacity (SWL) - 170 kg
- Max lifting capacity 300 kg
- Castors 1 x Ø125 mm, 2 x Ø150 mm
- Floor clearance 20/80 mm
- Table weight 240kg

ADJUSTMENT RANGES

- Hydraulic height adjustment (without mattresses) 670 - 1030 mm
- Hydraulic Trendelenburg / Reverse trendelenburg 28°/ 30°
- Hydraulic Lateral tilt (right/left) 20°/ 20°
- Gas spring assisted Back section (up/down) 65°/ 30°
- Gas spring assisted Leg section (up/down) 16°/ 85°, detachable
- Head section (up/down) 40°/ 40°, detachable

Made in Finland. CE.



WIDE RANGE OF OT-ACCESSORIES

