



















## MATERIAL AVAILABILITY

### Product-specific upholstery materials

	Capre FX treatment tables	Capre treatment tables	Delta treatment tables	Afrodite chairs	Special treatment tables (Bobath, Tilt, Multi-Tilt)	MTT benches	Capre BC chairs	Capre MC/RC chairs	Examination chairs 6800/6900	4000-series examination tables	Capre examination tables	Transportation stretchers	Field operating table	Lojer Antimicrobial-products	Working chairs
skai Fashion	•	•	•	•	•	•	•	•	•	•	•	•			•
skai Clinica	•	•	•	•	•		•	•	•	•	•	•			•
skai Conductive										•	•	•	•		•
Nature		•	•	•	•	•				•	•				•
Silvery														•	

## TECHNICAL SPECIFICATIONS

### Material/specification

																		
Material/ specification	Resistant to oil and fat	Antibacterial	Resistant to disinfectants (according to a separate listing)	Biocompatible	Waterproof	Resistant to secretions	Flame retardant	Resistant to abrasion	Lightfast	Tear-resistant	Long-lasting	Easy to clean	Electrically conductive	Antimicrobial	Fulfills Oeko-Tex 100 standard	Environmentally friendly	Fulfills REACH requirements	Fulfills toys norm DIN EN 71-3
skai Fashion			•				•	•	•	•	•	•					•	
skai Clinica	•	•	•	•		•	•	•	•	•	•	•					•	
skai Conductive			•						•	•	•	•	•				•	
Nature			•	•			•	•	•	•	•	•			•	•	•	•
Silvery	•	•	•		•	•	•	•	•	•	•	•		•		•	•	•

For more information: <http://www.lojer.com/service/material-options>