

# Alli Hydraulic Nursing Beds

An alternative to battery-operated nursing beds

**The Alli Hydraulic** is a hydraulically height-adjusted nursing bed with foot pedal operation. The Alli offers an alternative to battery-operated beds in environments where electricity use is limited, or where it poses a risk to patients. Other benefits of a hydraulic bed include reliability and lower costs over the life of the product. Apart from electric motors, the Alli Hydraulic features all that is expected of a modern nursing bed.

*"Durable, functional and mobile hydraulic nursing bed range."*

The base of the bed is made of hygienic ABS plastic with breathing holes. The foam mattress retains its qualities better on this kind of base than on a traditional mesh base. The base is easy to remove for cleaning and disinfecting.

Alli nursing beds allows stylish, quality interior décor without compromising functionality, coordinating perfectly with our range of bedside furniture and chairs. The bed ends are made of solid wood, giving a longer life than normal beds. The high-quality epoxy coating ensures a long service life; we are the only manufacturer to treat our steel parts with abrasive blasting before epoxy coating.



Alli H480 with side rails, mattress

**LOJER**® For easy care

Alli Hydraulic Nursing Bed	H-280	H-480	H-290	H-490
Number of sections	2	3	2	3
Width (incl. siderails)	90 cm	90 cm	100 cm	100 cm
Length 216 cm	216 cm	216 cm	216 cm	216 cm
Weight, kg (inc. siderails)	98	98	99	99
Height adjustment range, hydraulic, with foot pedals	41 .. 75 cm	41 .. 75 cm	41 .. 75 cm	41 .. 75 cm
Back section adjustment 0...70°, mechanical	X	X	X	X
Back section adjustment 0...70°, gas spring adjustment	•	•	•	•
Foot section adjustment, 12°, mechanical		X	-	X
Knee break adjustment 15°, mechanical		X	-	X
Knee break adjustment 15°, with gas spring		•	-	•
Lower body part adjustment, mechanical	•	-	•	-
Individually locking castors Ø125 mm	•	•	•	•
Centrally locking castors Ø125 mm	X	X	X	X
Wooden bed ends	X	X	X	X
Pair of double-locked wooden sides with locking button	•	•	•	•
Set of bumper rollers	•	•	•	•
Lifting pole with triangle	•	•	•	•
Lift support	•	•	•	•
IV-pole, straight or curved	•	•	•	•
Crutch stand	•	•	•	•
Padded coating for side rails	•	•	•	•
Urine bag holder	•	•	•	•
Reading and meal tray for placing on top of the bed sides	•	•	•	•
Recommended mattress size	205*78*13 cm	205*88*13 cm	205*78*13 cm	205*88*13 cm
Safe Working Load (SWL)	170 kg	170 kg	170 kg	170 kg
CE-marking	X	X	X	X
Warranty	2/10 years	2/10 years	2/10 years	2/10 years
Country of origin	Finland	Finland	Finland	Finland

Standard feature = X, Optional feature = •, Not available = -

## Accessories



Lifting pole and IV pole, curved



Lift support



Reading and meal tray for placing on top of the bed sides