

# Showering Furniture

Durable and Easy to Use

Lojer showering furniture has been developed in close co-operation with nursing staff, resulting in a range that is second to none.

The sturdy showering trolley features hydraulic height adjustment with foot pedals, with a range of 54 - 87 cm. The trolley is fitted with heat insulated padding and a pillow. The rails can be easily adjusted from both ends. The frame is epoxy-coated stainless steel. The trolley features high-quality durable, water-tight castors with 150mm wheels and a lockable direction wheel, making manoeuvring around the ward and corridors easy. The trolley tilts for drainage using a gas spring and has a drainage hose, complete with drip prevention and storage hook.

The 4080 shower chair is extremely sturdy and features hydraulic height adjustment, height adjustable foot supports which turn to the side, and lifting arm rests which can be used as a support frame. The frame of the chair is stainless steel and it features stainless, individually lockable 125 mm wheels and a locking direction wheel. The seat and back section are soft urethane with an anti-slip finish.

The hydraulically adjusted height range is 53 - 86 cm, and is operated by foot pedal. Fits over most standard toilet seats and can also be used in the sauna (with limitations, height up to 86 cm).

Shower Chair 4080



*"Developed in close co-operation with nursing staff."*



Shower Trolley 4310

**LOJER**® For easy care

Showering furniture	Shower trolley	Shower chair
Material, stainless steel, epoxy coated	X	X
Material, soft urethane with an anti-slip finish	-	X
Height adjustment, hydraulic, cm	54..87	53..86
Width, cm	66	48
Width max, cm	66 (77)	64
Weight, kg	74	47
Length, 190 cm	X	-
Depth, 45 cm	-	X
Back section length, 24 cm	-	X
Wheels, stainless steel, mm	150	125
Lockable direction wheel	X	X
Height adjustable foot supports which turn to the side	-	X
Lifting arm rests	-	X
Heat insulated pillow	X	-
Heat insulated mattress	X	-
Water pillow with hole	-	•
Toilet bowl and holder	-	•
Neck support	-	•
Lid for the seat	-	•
Plaster support kit	-	•
Safe Working Load (SWL)	190 kg	150 kg
CE-marking	X	X
Warranty 2/10 years	X	X
Country of Origin	Finland	Finland

Standard feature = X, Optional feature = •, Not available = -

