

Bobath Mat Table

A Wide Treatment Table for Neurological Therapy

Lojer Bobath mat treatment tables are designed primarily for neurological therapy and Bobath therapy. The tables are extremely sturdy and designed in such a way that the therapist can be on the table at the same time as the patient being treated. The table top is upholstered to a high standard and the sides are strengthened with plastic edging. A two-part table top with gas spring assisted, adjustable head section is available as an option.

"Classic and reliable model for professional use."

The construction of Lojer's Bobath table is extremely sturdy, and they require very little maintenance; all the moving joints, for example, feature maintenance-free sealed bearings. The mat table features electrically operated height adjustment between 44-88 cm and individually lockable 125 mm castors. The table offers high ground clearance, so a patient lift can fit under it easily.



**Bobath Mat Table with
two part table top (optional)**

LOJER® For easy care

Bobath Mat Table	
Table size, 120 x 200 cm	X
Table size 140 x 200 cm	•
Table size 105 x 200 cm	•
Weight, 92 kg (table size 120 x 200 cm)	X
Height adjustment, electric, range 44-88 cm	X
Foot control	X
Fixing plate for foot control	•
One part table top	X
Two-part table top (back section gas spring assisted)	•
Back section adjustment with gas spring, 0..75°	X
Breathing hole	•
Breathing hole padding (plug)	•
Individually locking castors, Ø125 mm	X
Centrally lockable castors, Ø125 mm	•
Set of supporting legs (lowest position, table standing on legs, not on castors)	•
Safe Working Load (SWL), 250 kg	X
Electric system, 230 V	X
SKAI Fashion upholstery (13 options)	X
CE-marking	X
Warranty 2 years	X
Country of origin	Finland

Standard feature = X, Optional feature = •



Bobath mat table with one-part table top