

Lojer Tilt Tables

Strengthening and Getting Patients Back on Their Feet

The tilt table is designed for helping patients get back on their feet and for strengthening therapy. The table features electrical height adjustment 45-80 cm, and electric angle adjustment from horizontal to 90 degrees.

Quality construction and features ensure patient safety. The table is ergonomically designed and the top is upholstered to a high standard, complete with cushioned edges. The hand grips can be locked into stainless steel rails at the sides of the table, and their position is steplessly adjustable. The product features two attachment straps as standard, as well as a wider belt to support the patient's lower limbs.

"Quality construction and carefully selected components ensure patient safety."



Lojer Tilt Table



Lojer Multi-Tilt Table

LOJER[®] For easy care

Features	Standard Tilt Table	Multi-Tilt Table
Length	190 cm	190 cm
Width	67 cm	67 cm
Weight	120 kg	128 kg
Angle adjustment range, lying surface (electric)	0...90°	0...90°
Angle adjustment range, back section (gas spring assisted)	-	0...80°
Angle adjustment speed, full movement (electric)	27 sec.	27 sec.
Safe working load (SWL), standard	200 kgs	200 kgs
Safe working load (SWL), battery operated	120 kgs	120 kgs
Universal side rails for accessories	X	X
Individually locking Ø 125mm castors	X	X
Centrally locking Ø 125mm castors	•	•
Choice of 13 SKAI upholstery colors	X	X
Manual emergency lowering for tilt function	•	•
Hand grips (2 ea)	X	X
Set of attachment belts (3 ea)	X	X
Battery operation (only with 230V)	•	•
Attachable work surface with padded elbow supports	•	X
Adjustable foot plate, one part	X	-
Adjustable foot plate, split	•	X
Hand grips with elbow support	•	•
Country of origin	Finland	Finland

Standard feature = X Optional feature = • not available = -

Accessories



Hand grips




Set of attachment belts



Working table



Foot plate, one part



Foot plate, split



Elbow support