

Capre F Treatment table

Versatility for physiotherapists

The modifiability, adjustability and silent operation of the sturdy Lojer Capre F improve functionality and enhance the comfort of use. Its multiple height settings, easy-to-use armrests and side supports, the stepless adjustable extension and flexion angle of the foot section improve the ergonomics and efficiency of the work being carried out. All Lojer Capre treatment tables come with high-quality, durable SKAI upholstery.

Lojer Capre – Leading the Way in Treatment Table Safety and Design



Lojer Capre F3



Lojer Capre F5



Scan the QR code to see the product video, or visit <http://youtu.be/TwpOU91j7IE>

LOJER[®] For easy care

Lojer Capre F	Capre F3	Capre F3H	Capre F5	Capre F5H
Table sections	3	3	5	5
Width (cm)	55,65,75	55,65,75	55,65,75	55,65,75
Upholstery SKAI Fashion (13 options)	•	•	•	•
Upholstery SKAI Clinica (3 options)	option	option	option	option
Anatomic design upholstery	option	option	option	option
Fully adjustable arm rests (gas springs)	-	-	•	•
Side supports with fast locking	option	option	option	option
Head support adjustment with gas spring	•	•	•	•
Double hinged head support	-	-	option	option
Leg raise with gas spring	+70°/-12°	+70°/-12°	+70°/-12°	+70°/-12°
Drainage with gas spring	option	option	option	option
Drainage with electric adjustment	option	-	option	-
Central locking castors (Ø75mm)	option	option	option	option
Foot adjustment (levelling)	•	•	•	•
Height adjustment range (cm)	45-92	49-90	45-92	49-90
Electric height adjustment with foot pedal	•	-	•	-
Hydraulic height adjustment with foot pedal	-	•	-	•
Additional hand control or foot pedal	option	-	option	-
Height adjustment with 360° foot bar	option	-	option	-
Paper roll holder (50 or 60 cm)	option	option	option	option
Battery operation	option	-	option	-
Safe working load (SWL)	210 kg	210 kg	210 kg	210 kg
CE-marking according to 93/42/EEC	•	•	•	•
Warranty (years)	2/10 y	2/10 y	2/10 y	2/10 y
Country of Manufacture	Finland	Finland	Finland	Finland

Explanation of markings: • =standard feature option =optional accessory - =not available



Capre F5 fully adjustable arm rests



Capre F3/F5 the stepless adjustable extension and flexion angle of the footsection



Drainage with gas spring or electric motor



Side supports with fast locking



Upholstery options F3, F5, F5 hourglass



See our upholstery options

<http://www.lojer.com/service/material-options>

LOJER® For easy care

Lojer Oy, Putajantie 42, P.O. Box 54, FI-38210 Sastamala, Finland
Tel. +358 (0)10 830 6700 | info@lojer.com | www.lojer.com