

Capre M Treatment Table

A design alternative for a basic treatment table

The sturdy, easy-to-use Lojer Capre M is a refreshing addition to the range of basic treatment tables in terms of its design and colour alternatives. Specifically intended for masseurs, the electrical Lojer Capre M has the best price-quality ratio on the market.

All Lojer Capre treatment tables have high-quality, durable SKAI upholstery. It is easy to keep the table clean, thanks to its straight surfaces and encased motor.

Lojer Capre – Leading the Way in Treatment Table Safety and Design



Lojer Capre M



Scan the QR code to see the product video, or visit <http://youtu.be/TwpOU91j7IE>

LOJER® For easy care

Lojer Capre M	Capre M2	Capre M2H	Capre M4	Capre M4H
Table sections	2	2	4	4
Width (cm)	55,60,65,80	55,60,65,80	55,60,65,80	55,60,65,80
Upholstery SKAI Fashion (13 options)	•	•	•	•
Upholstery SKAI Clinica (3 options)	option	option	option	option
Fully adjustable arm rests (gas springs)	-	-	•	•
Side supports with fast locking	option	option	option	option
Head support adjustment with gas spring	•	•	•	•
Double hinged head support	-	-	option	option
Central locking castors (Ø75mm)	option	option	option	option
Foot adjustment (levelling)	•	•	•	•
Height adjustment range (cm)	45-92	49-90	45-92	49-90
Electric height adjustment with foot pedal	•	-	•	-
Additional hand control or foot pedal	option	-	option	-
Height adjustment with 360° foot bar	option	-	option	-
Oval shaped frame tubing	•	•	•	•
Paper roll holder (50 or 60 cm)	option	option	option	option
Battery operation	option	-	option	-
Safe working load (SWL)	210 kg	210 kg	210 kg	210 kg
CE-marking according to 93/42/EEC	•	•	•	•
Warranty (years)	2/10 y	2/10 y	2/10 y	2/10 y
Country of Manufacture	Finland	Finland	Finland	Finland

Explanation of markings: • =standard feature option =optional accessory - =not available



Capre M4 fully adjustable arm rests



Side supports with fast locking



Paper roll holder



Central locking castors, 360° foot bar



Capre M2, Chili, 360° foot bar, castor set



Capre M4 Platin, 360° foot bar, castor set



Capre M4H, Ocean, hydraulic height adjustment



Capre M4, River, 360° foot bar, castor set, side supports



Upholstery options M2, M4



See our upholstery options

<http://www.lojer.com/service/material-options>