

# Capre F2 Treatment table

Sturdy Treatment Table for basic treatments and examination

The modifiability, adjustability and silent operation of the sturdy Lojer Capre F improve functionality and enhance the comfort of use. The Lojer Capre F2 features a two-piece table top. Table is designed for basic treatments and examination.

*High quality table for basic treatments and examination.*

All Lojer Capre treatment tables have a wide range of upholstery colors and several high-quality upholstery materials to choose from. It is easy to keep the table clean, thanks to its straight surfaces and encased motor.



**LOJER**<sup>®</sup> For easy care

Lojer Capre F2	
Table sections	2
Width (cm)	65, 70, 80
Length (cm)	196
Weight (kg)	95
Safe working load (kg)	210
Height adjustment range (cm)	45-92
Back section length (cm)	68
Back section adjustment	0-78°
Height adjustment with hand control	•
Electric back section adjustment with hand control	•
Battery use	option
360° height adjustment foot bar	option
Paper roll holder	option
Face hole cover	option
Foot pedal for height adjustment	option
Castors Ø75 mm	option
Electric system	24V
CE -marking	Yes
Warranty	2 years
Country of origin	Finland

Explanation of markings: • =standard feature option =optional accessory



Paper roll holder



Face hole cover



Hand control



See our upholstery options

<http://www.lojer.com/service/material-options>

**LOJER**® For easy care

Lojer Oy, Putajantie 42, P.O. Box 54, FI-38210 Sastamala, Finland  
Tel. +358 (0)10 830 6700 | [info@lojer.com](mailto:info@lojer.com) | [www.lojer.com](http://www.lojer.com)