

# ADA<sup>TM</sup>

## ADA bedside cabinet

Product pictured with  
optional accessories



The modular, durable and safe bedside cabinet Ada has been designed for heavy use in any hospital. With a wide range of accessories, the Ada cabinet can be modified to meet your needs.

The reading and dinner table is easy to adjust and can be quickly changed from side to side without any tools. It features a stepless height adjustment. Twin castors make the cabinet easy to move. Castor locks anchor the cabinet firmly in place. Lockable upper drawer helps to keep valuables safe. The frame is epoxy-coated steel, top and side panels as well as the dinner table is melamine. There's two drawers, of which the upper one is lockable. Top and side panels are grey and drawers white.

The angle of the reading and dinner table is adjustable. There's also a model without the reading and dinner table available. ADA bedside cabinet is CE marked.

**LOJER<sup>®</sup>** *For easy care*

Technical data and features	
Frame	Epoxy-coated steel
Height	90 cm
Width	47,5 cm
Height	90 cm
Depth	60 cm (113 cm with dinner table in front)
Table top dimensions	46 x 58,5 cm
Size of the dinner/reading table	38 x 48 cm (38 x 54,5 cm with holder)
Height adjustment of the dinner/reading table	780 - 1050 mm
Safe working load of the dinner/reading table	10 kg
Castors	Lockable twin 75 mm
Top and side panels	Grey melamine
Drawers	White, ABS-plastic
Table top cover (optional)	ABS-plastic
CE-marked	Yes

## Options and accessories



### Top and side panels

birch melamine ..... 10000 9187  
 beech melamine ..... 10000 9188

### Drawers

light blue fronts ..... 10000 9212



**Table top cover ABS**  
 plastic 10000 9262



**Slipper rack**  
 10000 9268



**Magazine basket**  
 10000 9266



**Bottle or waste bag rack**  
 10000 9288  
 patient intercom unit  
 10000 0629



**LOJER®** For easy care