

NOVA™

Lojer Merivaara NOVA bedside cabinet



The functional, safe and stylish Nova bedside cabinets are suitable for both hospitals and nursing homes. Thanks to its pull-through drawers, the cabinet can easily be used from both left and right-hand-side. Nova bedside cabinet is easy to handle and keep clean. Easy-rolling castors enable excellent and silent mobility. The color of the cabinet is light grey (melamine) and as an option also stained beech melamine is available.

Rounded edges and balanced construction ensure safe use also when using the reading / dinner table. The angle of the reading and dinner table is adjustable.

There's also a model without the reading and dinner table available.

Technical data and features	
Width	475 mm
Height	900 mm
Depth	600 mm, 1130 mm with dinner table in front
Height adjustment of the dinner/reading table	780 - 1050 mm
Size of the dinner/reading table	480 x 380 mm
Safe working load of the dinner/reading table	10 kg
Castors	Lockable twin 75 mm
Frame	Epoxy-coated steel
Dinner/reading table	ABS-plastic
Table top cover (optional)	White, ABS-plastic
Drawers	Grey melamine

Options and accessories



2660
stained beech melamine



2670
Model without the table



2680 (only the table)



Table top cover



LOJER® *For easy care*

© Lojer Group 2019. All Rights reserved | Lojer Group, Putajantie 42, 38210 Sastamala | Tel. + 358 (0) 10 830 6700 | info@lojer.com | www.lojer.com

The story of Lojer Group dates back to 1919 when Vammalan Konepaja was founded. In the early years, the company repaired f.ex. bicycles. Today, the company is the largest manufacturer of medical furniture in the Nordic countries and our products are used in 115 countries.