

Salli Fixed Height Hospital beds

Easy to use, modular and functional

Simple, easy to use, modular and functional Salli fixed height hospital beds were created for every day patient care. Made out of high quality raw materials and components these beds have 10 years warranty for welded parts and may stand very severe usage. These low-cost models gained popularity in several countries and became Lojer's bestseller for past years. Standard features include mechanical regulation of sections, removable bed end, ABS plastic surface. Optional side rails (2 types), lifting pole, IV stand, 125 mm. individual locking castors will add more to functionality of Salli fixed height beds.

"Easy to use, modular and functional bed collection."



Salli F280 Hospital Bed with castors

Salli F280 Hospital Bed with side rails, castors, mattress and IV-pole


LOJER® For easy care

■ Healthcare | Senior Care | Medical Wellness

SALLI	F180/190	F280/290	F380/390
Total width (incl. siderails)	90/100 cm	90/100 cm	90/100 cm
Total length	216 cm	216 cm	216 cm
Weight	70 kg	72 kg	74 kg
Height 54 cm (with castors 58 cm)	54/ 58 cm (without/with castors)	54/ 58 cm (without/with castors)	54/ 58 cm (without/with castors)
Safe Working Load (SWL)	200 kg	200 kg	200 kg
Back section adjustment 0...70° (mechanical)	-	X	X
Foot section adjustment 0/15° (mechanical)	-	-	X
Leg section adjustment 30° (mechanical)			X
Lower body part adjustment (mechanical)		•	
Gas spring assisted adjustment for back section	-	•	•
Gas spring assisted adjustment for knee-break	-	-	•
Individually locking castors Ø125	•	•	•
Metal bed end, curved/ straight, epoxy-coated/ chromium plated	X	X	X
Additional bed end curved/straight, epoxy-coated/chromium plated	•	•	•
Metal bed end, with tilting mechanism	•	•	•
Trendelenburg/Anti-trendelenburg 12°/14° (only for models with castors)	•	•	•
IV-pole straight model	•	•	•
IV-pole curved model	•	•	•
Lifting pole	•	•	•
Lift support	•	•	•
Pair of metal side rails (standard model) epoxy-coated or chromium plated, total length 178 cm	•	•	•
Pair of metal side rails, economy model (drop down model), epoxy-coated or chromium plated , total length 167 cm	•	•	•
Bumpers	•	•	•
Recommended mattress size 205 x 78/88 x 13 cm	•	•	•
CE-marking	X	X	X
Warranty	2/10 years	2/10 years	2/10 years
Country of origin	Finland	Finland	Finland

X = Standard - = not available • = Optional


Accessories



IV pole



Side rails



Lifting pole