

Capre RC Recovery Chair

For recovery after outpatient surgery



Beautifully designed patient chair for recovery after outpatient surgery – also suitable for dialysis and cytostatic treatments.

The luxurious Capre multipurpose chair has been designed together with surgeons to fulfil the demands of today's healthcare processes in post-day surgery recovery. The chair is used instead of the traditional hospital bed in order to improve patient vitality during recovery. The Capre chair features pleasantly contoured, soft imitation leather, individual adjustment possibilities (with five motors) and built-in comforts, such as a tablet pc stand and reading light. All of this combines to improve patients' quality of life during recovery, whilst also improving the quality image of the healthcare facility.

The control system allows the independent adjustment of the chair height, back section, seat angle and leg sections. The leg section is split and motorised. Both the head section and the split leg section also feature a length adjustment mechanism. The chair features five pre-programmed positions: three treatment positions (sitting, beach chair and lying), Trendelenburg position and one memory slot for storing the customer's own settings. Due to its modular design, the chair is also suitable for dialysis and cytostatic treatments, and as an examination chair.



Scan the QR code
to see the product video, or visit
<http://youtu.be/EbrNPcGUhyk>

LOJER® For easy care

Technical specifications

- **Total length:** 186..213 cm
- **Height adjustment:** 60..85 cm
- **Castors:** Ø 75mm, centrally locking
- **Back section adjustment:** 0°..+90°
- **Head section adjustment:** 0°..10°
- **Back section length:** 600 mm
- **Seat adjustment:** +12°..0°
- **Seat width:** 620 mm
- **Total width:** 880 mm
- **Leg section adjustment:** -80°...0°
- **Arm support adjustment:** 0°..180° rotation, -45°..+20° angle
- **Safe Working Load (SWL):** 160 kg
- **USB output:** Power 4W
- **Electric system:** 24 V
- **Warranty:** 2 years
- **Country of Manufacture:** Finland



Features	
CE -marked medical device	x
Luxurious anatomically contoured upholstery	x
Two colour options for the plastic parts	x
Upholstery SKAI Fashion (16 options)	x
Design castors for moving the chair, centrally lockable	x
Intelligent control system with pre-programmed positions	x
Trendelenburg (12° max.)	x
Motorised adjustment of back section	x
Motorised adjustment of seat angle	x
Motorised height adjustment	x
Motorised split leg section with length adjustment	x
Head section length and angle adjustment	x
Arm rests stay horizontal regardless of back section angle	x
180° rotation of arm rests, angle and height adjustment	x
Safe, quiet 24V electrical system, USB charger	x
Warranty 2 years	x
Optional features:	
Coffee tray	o
Tablet PC stand with flexible neck	o
Face cut-out	o
Plug for face cut-out	o
Decorative stitching in the upholstery	o
LED light (usb)	o
Infusion stand	o
Foot controller for hi-low and back section	o
Foot pedal for Trendelenburg	o
Wired hand control (instead of flexible arm)	o

Explanation of markings: x = standard feature, o = optional accessory, - = not available

Accessories



Coffee tray



Tablet support
with flexible arm



Stitching



Led light



Foot controller