

Capre E1 and E2 Examination Tables

Modern and unique look for demanding projects



Product pictured with optional accessories.



Main features:

1. Quiet and easy to use
2. Stands firmly on rubber paws
3. Easy to keep clean


The Lojer Capre E1 and E2 examination tables are in a class of its own when it comes to design and function. There are two Lojer Capre Examination table models. The Lojer Capre E1 examination table features a one-piece tabletop and the Lojer E2 examination table a two-piece tabletop. Both are suitable for general examination use. Low minimum height and quiet motor add to their appeal in small examination room environments.

Unique design will delight also demanding customers and will surely upgrade the interior of all kind of medical premises. A large number of examination tables are typically placed in examination rooms with no need for manoeuvring – large castors are therefore unnecessary. Moving the table during cleaning is possible with the help of small castors, which are available as an accessory. When table is in use, it stands firmly on its own rubber feet and is therefore extremely sturdy.

Table's great height adjustment range allows the best possible ergonomics for healthcare personnel, not to mention for the patient's high comfort. With the lowest height of 45 cm it is easy for older people and for patients on the wheelchair to get on the table. Height adjustment is convenient with an optional 360° height adjustment foot bar.

LOJER® For easy care

Technical data and features	E1	E2
Total width	69 cm (optional 84 cm)	68 cm (optional 83 cm)
Patient surface width	65 cm (optional 80 cm)	65 cm (optional 80 cm)
Total length	196 cm	196 cm
Height adjustment	45 - 92 cm	45 - 92 cm
Back section length	-	72 cm
Seat and leg section length	-	123 cm
Back section adjustment	-	0 - 78°
Total weight	95 kg	95 kg
Safe working Load (SWL)	210 kg	210 kg
Castors (optional)	Ø75 mm	Ø75 mm
Upholstery	Several color and material options, please ask for more info	Several color and material options, please ask for more info
Hand control	Optional	Standard
Foot control	Standard	Optional
Electric system	24V	24V
CE-marking	Yes	Yes
Warranty	2 years	2 years
Country of manufacture	Finland	Finland



Factory installed options:

- Foot control for height adjustment (E2)
- Hand control for height adjustment (E1)
- 360° height adjustment foot bar (incl. safety switch)
- Safety switch
- Patient surface width 80 cm
- Central locking castors, ø 75mm
- Mounting for IV-pole (pair)

Other accessories:

- IV-pole
- Side supports for E2
- Paper roll holder and cutter
- Neck cushion



Capre E1 Examination Table



LOJER® For easy care