## Modux Folding Care Bed

Best-selling folding care bed



## Main features:

- 1. Small footprint for storage (just 90/100 x 60 x 127 cm when folded for transport or storage)
- 2. Folds using electric motors at the press of a button, no tools required
- 3. Available in 80 cm and 90 cm widths

The patented Modux care bed folds quickly and easily for transport or storage. The Modux is the perfect care bed for units loaning equipment for home care, or as a flexible solution for nursing homes. Although the Modux is a folding bed it still has all the same features as other nursing beds, including a very low minimum height (from 28 cm). Patented Modux nursing beds are flexible and versatile, providing a pleasant care experience for the patient and ergonomic working for care staff. The Modux can also be adjusted to a comfortable position similar to a beach chair.

**LOJER**<sup>®</sup> For easy care

Senior Care | Healthcare | Medical Wellness




Technical data	
Transportation dimensions (w x l)	90/100 x 60 x 127 cm
Overall width with rails	90 or 100 cm
Overall length	223 cm
Height adjustment range	30 - 69 cm (125 mm castors) 28 - 67 cm (100 mm castors)
Recommended mattress size	78-80 or 88-90 x 205 x 9-13 cm
Total weight	92 or 96 kg
Castors	Ø125 mm locking castors
Safe Working Load (SWL):	170 kg
Electrical system	24V
Back section adjustment	0-71°
Thigh section adjustment	0-38°
Calf section adjustment	0-12°
Anti-trendelenburg	0-12°
Other standard features	Hand control   solid birch bed-ends (pair)   double-locking wooden siderails with locking button (pair)   integrated fixing for lifting pole (head-end)   integrated fixing for IV-pole (head-end)
CE -marked	Yes
Warranty	2 years
Country of manufacture	Finland

## Factory installed options:

• Ø100 mm castors

## Other accessories:

- Lifting pole with triangle handle
- Lift support
- IV-pole with one hand quick adjustment
- Reading and dining tray, on top of the side rails
- Reading lamp
- Padded coating on top of the side rail, one unit
- Bed ends: birch







© Lojer Group 2020. All Rights reserved | Lojer Group, Putajantie 42, 38210 Sastamala, Finland | Tel. + 358 (0) 10 830 6700 | info@lojer.com | www.lojer.com The story of Lojer Group dates back to 1919 when Vammalan Konepaja was founded. In the early years, the company repaired f.ex. bicycles. Today, the company is the largest manufacturer of medical furniture in the Nordic countries and our products are used in 115 countries.