

ScanAfia XL nursing bed

A Quality New Arrival for Home and Home-like Care



The large twin castors make the bed easy to manoeuvre yet extremely sturdy when locked in position.

The special sized ScanAfia XL nursing bed has been designed with larger and heavier patients in mind. The bed is particularly suitable for nursing homes, where the average size of patients has increased. The special sized bed is 120 cm wide, with a maximum patient weight of 250 kg.

The ScanAfia nursing bed can be combined with our range of tables and chairs, creating a harmonious décor without compromising functionality. The solid white oak bed ends give the bed a longer life cycle than standard beds. Our high quality epoxy coating also contributes to the long lifespan of our beds – we are the only manufacturer to abrasive blast steel components before painting.

The bed features a safe 24v low voltage electrical system, and its electrical motors and control units are splash resistant. The beds can be fitted with a separate nurse panel as an optional extra, allowing the functions of the bed to be restricted to nurse control only if necessary. The design of the bed is based on the 60601-1-52 nursing bed standard and the range has been impartially tested in an independent laboratory.

LOJER[®] For easy care

TECHNICAL DATA

Overall width with rails	130 cm
Overall length	216 cm
Wooden parts	Solid wood bed ends (white-oak), Two wooden sides with locking button
Lying surface	Bed section base IKI high pressure laminate
Steel parts	Steel parts abrasive blasted and epoxy coated
Height adjustment range	32..80 or 34..82 cm
Recommended mattress size	205x 118 x 13 cm
Total weight (inc. rails)	145 kg
Castors	Ø150 mm centrally lockable double castors and direction wheel
Safe Working Load (SWL)	285 kg
Electrical system	24 V, Electrical parts splash-proof and double insulated for safety
Back section adjustment	70°
Thigh section adjustment	12°
Calf section adjustment	15°
CE-marked	Yes
Warranty	General warranty 2 years, steel frame 10 years
Country of Manufacture	Finland



See our options

<http://www.lojer.com/service/material-options>

LOJER® *For easy care*

Lojer Oy, Putajantie 42, P.O. Box 54, FI-38210 Sastamala, Finland
Tel. +358 (0)10 830 6700 | info@lojer.com | www.lojer.com