

ScanAfia XS Hospital Bed

Extremely low and easy to move

Product pictured with optional accessories.



Main features:

1. Ergonomic working position thanks to impressive height adjustment range
2. High safe working load (SWL): 260 kg
3. Easy to manoeuvre thanks to its double castors (optional)

The quality ScanAfia XS hospital bed is the safe choice for hospital wards. The bed has been designed for use by healthcare professionals in hospitals, healthcare centres and similar facilities, and it fulfils the requirements of the new EN 60601-2-52:2010 standard for medical beds.

Although it can take loads up to 260 kg (SWL), the ScanAfia XS hospital bed is still light and easy to manoeuvre, thanks to the optional choice of double castors. The design places emphasis on both ergonomic working and patient comfort and accessibility; its exceptionally low minimum height makes work easier for nurses and makes getting out of bed easier for patients. The back section slide mechanism also improves patient comfort. The bed also features an extremely broad Trendelenburg / Anti-Trendelenburg adjustment range.

LOJER® For easy care


Technical data and features	ScanAfia X ICU
Total width (incl. siderails)	89,5 or 99,5 cm
Total length	215 cm
Height adjustment	34 - 82 cm (125 mm castors)
Recommended mattress size	32 - 80 cm (125 mm twin castors)
Total weight (without rails)	34 - 82 cm (150 mm twin castors)
Castors	205 x 78/88 x 13 cm
Safe Working Load (SWL)	80 cm wide 108 kg, 90 cm wide 110 kg
Electrical system	Ø125 centrally lockable wheels (Ø125 mm and Ø150 mm twin castors as an option)
Trendelenburg/Anti-Trendelenburg	260 kg
Back section adjustment	24V
Thigh section adjustment	12°/ 20°
Calf section adjustment	70°
Other standard features	Hand controller manual emergency lowering of back section back section slide lowest height starting from 32 cm vertical upward moment reduces the need for space removable bed ends, available in standard colours white, blue, yellow and light green chrome plated or painted side rails; normal and long versions fixing for restraining belt bumper rollers built-in place for lifting pole and two IV poles
CE-marked	Yes
Warranty	2 years
Country of Manufacture	Finland

Factory installed options:

- Battery driven
- 5th wheel to assist moving
- 125 mm twin castors
- 150 mm twin castors

Other accessories:

- Special length (extendable +15 cm) model
- Extra height side rails
- Lifting pole with triangle handle
- Lift support
- Lift support version 2
- IV pole with one hand quick adjustment
- Foldable bed end (works also as an auxiliary tray for linen)
- Nurse locking panel (allows functions to be locked to prevent patient use)
- Electric CPR
- Urine bag holder
- Reading and dining tray, on top of the side rails
- Bed end (for head end of the bed)
- Monitor tray / end
- Flexible hand controller holder



LOJER® For easy care