

Bedside Cabinet 2300

Lojer Bedside Cabinet 2300 is a trendy novelty for home-like care. Available in two stylish finishes: white oak and beech. The cabinet includes a drawer with optional locking, and single-locking castors are also available as accessories.



Bedside Cabinet:
white oak



Technical data and features	
Dimensions table top (w x d)	40 x 50 cm
Dimensions table frame (w x d)	38 x 46.5 cm
Height	65 cm (72.5 cm with castors)
Weight	~19 kg
Drawer dimensions (inside)	37 x 29 x 12 cm
Warranty	2 years
Country of manufacture	Finland



Features:

- Color options include white oak and beech
- One drawer and a shelf

Accessories:

- Lockable castors
- Lock for the drawer
- Beech or white-oak door