



FYSICON EVAC MiniMAX

Antidecubitus evacuation

The new anti-decubitus mattress with integrated evacuation bands

Purpose of use

- Fysicon Evacuation MiniMAX is a pressure-reducing mattress for patients with low / middle risk of pressure ulcers in hospitals, nursing homes and home care.
- MiniMAX can be used for patients 40-150 kg and it is giving the best support and softness to all users. If using by a heavy person or 24 hours/day, lifetime is shorter than normally (recommendation in public health care 3-5 years).

Product structure

- MiniMAX mattress core has 3 layers: on top the soft viscoelastic foam, in the middle and bottom more supporting and heavier layers.
- Surface and sides of the cover is polyurethane coated hygienic fabric, which is impermeable to liquids but still breathable. Material is bi-elastic and treated against bacteria and fungi (active biocidal substance is zinc pyrithione (EU 528/2012). Zipper on two sides, protected with flap. The bottom is made of durable and slippery material with integrated evacuation bands: there are 2 pcs of bands for patient fixing and pulling bands in the both ends.

Benefits

- The same mattress is suitable for patients with different sizes because the surface is soft and bottom layers are supporting heavier users.
- Patient's condition and care are easier when bedsores can be avoided and change of position is not needed so often.
- Because of the integrated evacuation bands the same mattress is easy to use in the normal patient care and in the rescue situation. Simple to use and easy to bed when there is not any need to a separate rescue sheet under the mattress.



Cleaning and maintenance

- Disinfection by wiping according to the instructions with diluted chlorine-amine or an alcohol-based detergent or ethanol /chlorine max 1000 ppm, rinse and dry thoroughly. If necessary, wash the cover 60°C inside out, gentle tumble drying at low temperature.
- Turn the mattress regularly to ensure even wear. Replace broken cover to a new one.
- The bed frame should be lathed or perforated to allow for proper ventilation.

Flame retardancy & CE

Flammability by standard IMO FTPC part 9. Fysicon mattresses are CE-marked. Country of origin Finland.

Size

Standard sizes 75*200*13, 80*200*13, 85*200*13 and 90*200*13 cm.

Other sizes on request.