

Capre Examination Tables

Modern and unique look for demanding projects



Product pictured with optional accessories.

Main features:

1. Modern and unique looking design for demanding projects
2. Stands firmly on non-marking rubber feet, making the table extremely sturdy
3. Good number of factory-installed options and accessories available


The Lojer Capre Examination tables are in a class of their own when it comes to design and function.

The unique design upgrades an interior of all kinds of medical premises. Examination tables are often placed in examination rooms with rarely a need to maneuver. With the Capre EX models it is possible to choose optional 75 mm castors that improve movability for example for cleaning purposes.

The height adjustment gives best possible ergonomic for professionals and for the patients' comfort. The lowest height of 45 cm is easy for older patients and patients using a wheelchair. Capre EX – examination tables have a range of accessories from side supports to anesthesia frame and a sampling arm rest. Optional 360°- frame control that gives convenience for height adjustment function improving the user experience.

This model is available as Capre EX1 that has a one-part table-top. Capre EX2 model with a two-part table-top and the hydraulic model Capre EX2H with a two-part table-top.

Technical data	EX1	EX2	EX2H
Total width	70 cm (80 cm optional)	70 cm (80 cm optional)	70 cm (80 cm optional)
Patient surface width	70 cm (80 cm optional)	70 cm (80 cm optional)	70 cm (80 cm optional)
Total length	200 cm	200 cm	200 cm
Patient surface length	200 cm	200 cm	200 cm
Total weight	100 kg	100 kg	95 kg
Height adjustment	45 – 92 cm, electric	45 – 92 cm, electric	49 – 90 cm, hydraulic
Back section length	-	70 cm	70 cm
Back section adjustment	-	0-73°	0-73°
Central locking castors (optional)	ø 75 mm	ø 75 mm	ø 75 mm
Other standard features	Electric height adjustment with hand control Rubber feet One adapter for IV-pole and anesthesia frame (6 fixings for installation)	Electric height and back section adjustment with hand control Rubber feet One adapter for IV-pole and anesthesia frame (6 fixings for installation)	Gas spring adjustable back section adjustment Hydraulic height adjustment with foot pedals Rubber feet One adapter for IV-pole and anesthesia frame (6 fixings for installation)
Upholstery	Several colors and material options, see the Upholstery card or lojer.com for more info	Several colors and material options, see the Upholstery card or lojer.com for more info	Several colors and material options, see the Upholstery card or lojer.com for more info
Safe Working Load (SWL)	180 kg	210 kg	180 kg
CE marked	Yes	Yes	Yes
Warranty	2 years	2 years	2 years
Country of Manufacture	Finland	Finland	Finland



Factory installed options:

- Foot control for height adjustment to EX1 (incl. hand control)
- Foot control for height and back section adjustment to EX2 (incl. hand control)
- Gas spring adjustable back section adjustment to EX2 (incl. hand control for height adjustment)
- 360° height adjustment Foot control for EX1 and EX2 (incl. safety switch and hand control)
- Safety switch
- Patient surface width 80 cm
- Central locking castors, ø 75 mm

Other accessories:

- IV-pole
- Anesthesia frame
- Extra adapter for IV-pole and Anesthesia frame (6 fixings for installation)
- Accessory rail (max. 6 pcs)
- Examination light with Lojer adapter (6 fixings for installation)
- Examination light with rail clamp (can be attached to accessory rails on the sides of the table)
- Neck cushion
- Paper roll holder
- Sampling arm rest
- Side rails

