

Lojer Delta treatment tables

Classic, hard-wearing, user-friendly favourites



LOJER®

For Easy
Care





LDP5



LDP4



LDP3



LDP2

Lojer Delta Professional

The Lojer Delta Professional treatment tables are designed with the needs of healthcare professionals in mind, offering functionality, durability, and safety. The extensive product development behind these tables ensures they remain stable and sturdy, regardless of working height or position.

Key features:

- Electric height adjustment from all sides
- Double hinged (option) head section with gas spring adjustment
- Armrests with user-friendly height and angle adjustment
- Manual or electric drainage positioning and footrest angle adjustment to 68°
- Standard width 55 cm with customisation to 50 cm or 60 cm available
- Safe working load 210 kg



Tables come with a variety of accessories and a wide choice of options from our Future and Nature upholstery range (link in QR code).



LDS5



LDS4

Lojer Delta Standard

Delta Standard is our series of basic treatment tables that are more affordable than the Professional series. All the basic features that a quality-conscious user wants can also be found in our Standard series.

Key features:

- Electric height adjustment from all sides
- Gas-spring adjustable head section
- Armrests with user-friendly height and angle adjustment
- Manual drainage positioning and footrest angle adjustment to 68°
- Safe working load 210 kg
- Hydraulic models also available



Tables come with a variety of accessories and a wide choice of options from our Future and Nature upholstery range (link in QR code).

Technical data and features

Delta Professional	LDP2	LDP3	LDP4	LDP5
Width	50, 55, 60 cm	50, 55, 60 cm	50, 55, 60 cm	50, 55, 60 cm
Total length	205 cm	205 cm	205 cm	205 cm
Weight	98 kg	98 kg	98 kg	98 kg
Height adjustment range	54-99 cm (with extended height adjustment range 45-99 cm)	56-102 cm (with extended height adjustment range 48-102 cm)	54-99 cm (with extended height adjustment range 45-99 cm)	56-102 cm (with extended height adjustment range 48-102 cm)
Leg section adjustment	–	0-68°	–	0-68°
Other features	Foot control/frame control/hand control O- or U-shaped head section	Manual/electrical drainage Manual leg section adjustment Foot control/frame control/hand control O- or U-shaped head section	Foot control/frame control/hand control O- or U-shaped head section Arm rests	Manual/electrical drainage Manual leg section adjustment Foot control/frame control/hand control O- or U-shaped head section Arm rests
Electric system	220-240V	220-240V	220-240V	220-240V
Delta Standard	LDS4	LDS4H	LDS5	LDS5H
Width	50, 55, 60 cm	50, 55, 60 cm	50, 55, 60 cm	50, 55, 60 cm
Total length	205 cm	205 cm	205 cm	205 cm
Weight	75 kg	59 kg	75 kg	59 kg
Height adjustment range	54-85 cm	54-85 cm	57-88 cm	57-88 cm
Leg section adjustment	–	–	0-68°	0-68°
Other features	Foot control/frame control/hand control	Hydraulic height adjustment, foot pedal O- or U-shaped head section Arm rests	Manual drainage Manual leg section adjustment Foot control/frame control/hand control O- or U-shaped head section Arm rests	Hydraulic height adjustment, foot pedal Manual drainage Manual leg section adjustment O- or U-shaped head section Arm rests
Electric system	220-240V	–	220-240V	–
Professional and Standard, data for all models				
Upholstery	Several colours and material options, see the Upholstery card for more info			
CE Marked	Yes (The device is a CE-marked medical device according to MDR 2017/745 class 1)			
Safe Working Load (SWL)	210 kg			
Warranty	2 years			
Country of Origin	Finland			
Manufacturer	Lojer Oy			

Accessories



Wheel mechanism



Side supports



Height adjustment bar



Paper roll holder



Electrical drainage



Double-hinged face section

