

CARENA PLUS

Quality hospital bed



Product pictured with optional accessories.



Options:



Main features:

1. Clear and improved structure
2. Nurse Control Panel as an option
3. Mechanical and electrical Trendelenburg adjustments
4. Safe and easily removable side rails

Carena Plus hospital bed has a clear and improved structure enabling fast cleaning and improved hygiene management of the frame. Designed with unparalleled ergonomics for patients and staff. The improved height range from the lowest 34 cm position enables easy access for smaller patients. Yet, this hospital bed has high safe working load (SWL) of 330 kg and is still very light to manoeuvre.

Carena Plus fulfills all requirements for medical devices: IEC 60601-2-52, IEC 60601-1-2 (EMC) and SFS-EN 60601-1

Optional Nurse Control Panel enhances the usability of the bed ensuring various positions and memorized positions. Rapid locks for locking the bed ends, as a standard feature, ensure easy and fast operation in emergency situations.

This hospital bed model is available as a washable model with a wash cycle counter.

Technical data and features	
Total width (with siderails)	90 / 100 cm
Patient surface width	80 / 90 cm
Total length	215 cm (extended 245 cm)
Patient surface length	202 cm (extended 232 cm)
Height adjustment	34–84 cm (Ø150 mm castors)
Recommended mattress size	75 / 85 x 200 x 13 cm
Total weight (without rails)	80 cm 90 kg / 90 cm 110 kg
Castors	Ø150 central locking castors: 2 pedals and a brake pedal at foot end
Safe Working Load (SWL)	330 kg
Electrical system	IPX6
Back section adjustment	0–70°
Thigh section adjustment	0–29°
Calf section adjustment	0–24°
Trendelenburg / anti-Trendelenburg	12° / 6° electric / manual
Removable bed ends	White, blue, beech, birch, green
Side rails	Different options available
Other features	Hand control, frame color; painted silver; mattress base cover ABS-plastic (ventilated and detachable) bumpers, fixing for restraining straps, vertical raising motion, simultaneous movements of different parts (f.ex.height and back section)
CE marked	Yes
Warranty	2 years
Country of Manufacture	Finland (Lojer Oy)

Factory installed options:

- CPR release handle
- Extendable +30 cm leg section (for 80 and 90 cm patient surfaces), ABS-plastic
- Adjustable metal net cover
- Mattress base for X-ray
- 5th wheel
- Different castor sizes
- Battery
(comes always with a nurse control panel)
- Potential equalisation system
- Extra hand control

Other accessories:

- Nurse control panel
- Bed end detachable
Frame tube: chromed and painted
- Bed end detachable and foldable.
Frame tube: chromed and painted
- Accessory basket 1-3: small, large and bed end options
- Infusion rod with 4 hooks or attachment to the lifting pole
- Drainage bag holder
- Oxygen bottle holder 5 kg
- Accessory rails
- Hook for electrical cable
- Grab handle
- Lifting pole
- Adapter for grab handle

