

# Shower Trolley 4310

Durable and easy to use



The sturdy showering trolley features hydraulic height adjustment with foot pedals, with a range of 54 - 89 cm. The trolley is fitted with heat insulated padding and a pillow.

The frame is epoxy-coated stainless steel. The trolley features high-quality durable, water-tight castors with 150 mm twin castors and a lockable direction wheel, making manoeuvring around the ward and corridors easy. The trolley tilts for drainage using a gas spring and has a drainage hose, complete with drip prevention and storage hook.

## Technical specifications and features

Main material	Stainless steel, epoxy coated
Height adjustment	54–89 cm
Patient surface width	66 cm
Total width	78 cm
Patient surface length	190 cm
Total length	204 cm
Subframe width	72 cm
Subframe length	114 cm
Rail adjustment angle	0–190°
Drain position angle	0–6°
Weight	74 kg
Safe Working Load (SWL)	190 kg
CE marked	Yes (The device is a CE marked medical device according to MDR 2017/745 class 1)
Warranty	2 years
Country of manufacture / Manufacturer	Finland / Lojer Oy

