

FUTURA PLUS™

Lojer Merivaara Futura Plus hydraulic hospital bed



Product pictured with optional accessories



Main features:

1. Sliding and ergonomically supportive back section
2. User-friendly hydraulic adjustment
3. Comfortable and reliable basic hospital bed

A true bestseller, the classic Futura Plus bed combines comfort, intuitive ease of use, and long-lasting reliability. It is a durable patient bed for hospitals, senior care units and nursing homes.

The Futura Plus complies with the medical bed standard IEC 60601-2-52. The Futura Plus bed is available in hydraulically adjustable models. The bed is easy and comfortable to use. Sliding and ergonomic back section provides better comfort for the patient.

The range:

- 2-part bed section (200 x 81 cm), with hydraulic height and back section adjustment
- 3-part bed section (200 x 81 cm), with hydraulic height, back and leg section adjustment



As a safety feature, the bed has gas-spring assisted functions for Trendelenburg adjustment.



Merivaara's Futura Plus has been appreciated by users in several countries for over 20 years.



There are also optional accessories which help the patient move in and out of the bed by him/herself.



The accessories attached to the head end do not move with the back section.



The new sliding back section provides better comfort for the patient. The bed comes with removable and easy-to-clean mattress base plates made of ABS plastic.



The back section can be adjusted with a gas-spring assisted handle.

LOJER®

For Easy
Care

Technical data and features	
Total length	2225 mm
Total width	912 mm
Patient surface length	812 mm
Patient surface width	2001 mm
Height adjustment	450–840 mm
Recommended mattress size	80 x 200 cm
Castors	Ø150 mm central locking castors: 2 pedals and a braking bar at foot end
Maximum patient weight	210 kg
Safe Working Load (SWL)	230 kg
Trendelenburg/Anti-Trendelenburg	26° / 12°
Back section adjustment	0–70°
Leg section adjustment	0–38°
Other features	Mattress base cover ABS-plastic (ventilated ja detachable) bumpers fixing for restraining straps vertical raising motion simultaneous movements of different parts (f.ex. height and back section)
CE marked	Yes (The device is a CE marked medical device according to MDR 2017/745 class 1)
Warranty	2 years
Country of Manufacture	Finland (Lojer Oy)

Factory installed options:

- Bed end, detachable. Frame tube: chromed or painted.
Bed end plate color white. Other colors: light blue, light green, beech and birch
- Bed end, detachable and foldable (works as auxiliary tray for linen).
Frame tube: chromed or painted. Bed end plate color white.
Other
colors: light blue, light green, beech and birch
- Side rails, chromed 1-sectional foldable
- 3/4 side rails, chromed foldable
- Adapter for grab handle
- Trendelenburg
- 4 braking pedals
- Central locking Ø150mm twin castors
- Lifting pole
- Grab handle
- Oxygen bottle holder 5kg

Other accessories:

- Journal basket
- Monitor tray
- Infusion rod with 4 hooks or attachment to the lifting pole
- Drainage bag holder
- Urine bottle holder

