

NOVA™

Lojer Merivaara NOVA bedside cabinet



The functional, safe and stylish Nova bedside cabinets are suitable for both hospitals and nursing homes. Thanks to its pull-through drawers, the cabinet can easily be used from both left and right-hand-side. Nova bedside cabinet is easy to handle and keep clean. Easy-rolling castors enable excellent and silent mobility. The color of the cabinet is light grey (melamine) and as an option also stained beech melamine is available.

Rounded edges and balanced construction ensure safe use also when using the reading / dinner table. The angle of the reading and dinner table is adjustable.

There's also a model without the reading and dinner table available (model 2670).

Technical data and features	
Frame	Anodized aluminium profile and composite material
Total width (dinner table turned down)	66 cm
Table top height	86,2 cm model 2670: 79,4 cm
Table top width	48 cm
Drawer (w x d)	46 x 39 cm
Table top plate (w x d)	48 x 37,5 cm
Dinner/reading table	36 x 53 cm
Height adjustment of the dinner/reading table	78–106 cm
Castors	Individually locking twin 75 mm model 2670: Individually locking twin 50mm
Safe working load (SWL) of the dinner/reading table	10 kg
Weight	33 kg model 2670: 25 kg
Dinner/reading table	Plastic composite
Drawers	Melamine
Table top cover (optional)	ABS-plastic

Options and accessories



2660
stained beech melamine



2670
Model without the table



2680 (only the table)



Table top cover

The top drawer is also available in a lockable version.

