

# 4040X General examination table

Safe and versatile examination table



Product pictured with optional accessories.

## Main features:

1. Exceptionally wide height adjustment range 40–95 cm
2. An optional two-part hands-free foot control (separate adjustment for height and backrest), together with the table's clear structure, ensures a high level of hygiene
3. A new safety feature detects any pinching situation, stops the movement, and raises the table to release the pressure

The 4040X general examination table is suitable for use in various healthcare facilities. Thanks to the table's design and exceptionally low minimum height, both children and elderly patients can easily get onto the table. An optional lift support further enhances safety and comfort during access. The table is lightweight and easy to move. It features advanced safety functions: it automatically stops and rises when movement is obstructed, it can be locked using any of the controls, and a manual emergency backrest lowering function is available as an option.

The exceptionally wide height adjustment range and high load capacity enable ergonomic working positions and enhance safety during patient examinations and transfers. With the optional two-part hands-free foot control, the height and backrest of the table can be adjusted separately from either side using the foot pedal, freeing the professional's hands for patient care. The new digital control unit of the examination table notifies when maintenance is required, and with remote support, customers can quickly receive assistance when needed.

Also available as the hydraulic 4040XH version, where the height adjustment is operated by a foot pedal, and the backrest and Trendelenburg positions are adjusted with gas spring assistance.

Technical data and features	
Total width	75 cm (80/90 cm optional), 4040XH model: optional 80 cm
Patient surface width	70 cm (80/90 cm optional), 4040XH model: optional 80 cm
Total length	200 cm
Patient surface length	200 cm
Total weight	108 kg
Height adjustment range (depends on castor size and product model)	40–95 cm (Ø125 mm castors)
Back section length	70 cm
Back section adjustment	73°
Electric Trendelenburg / anti-Trendelenburg	20° (optional)
Manual Trendelenburg	12° (optional)
Castors	Ø125 mm (Individually lockable)
Castors (optional)	Centrally lockable Ø125/150 twin castors
Other features	Maximum of 6 different mounting positions for the fluid tray   Table can be locked from all controls   Safety feature: the table automatically stops and rises when movement is obstructed
Upholstery	5 cm or, optionally, thicker 7 cm upholstery (*)
Safe working load (SWL)	225 kg
CE marking	Yes (This device is a CE-marked medical device, MDR 2017/745 class 1)
Warranty	2 years
Country of Manufacture	Finland (Lojer Oy)

#### Factory installed options:

- Hand control
- A hands-free adjustment bar for height and back section
- Single- or two-part foot control
- Manual Trendelenburg
- Electric Trendelenburg and anti-Trendelenburg, including battery
- Battery (emergency use)
- Centrally lockable 125/150 mm twin castors
- Facehole
- Manual emergency lowering of the backrest

#### Other accessories:

- Side rail
- IV-pole
- Arm rest
- Neck cushion
- Examination light
- Push handles
- Adapter for accessories
- Accessory rail  
380 mm x 10mm x 25 mm (max. 6 units)
- All accessories fitted to accessory rail fixing
- Anaesthesia frame
- Lift support
- Collision bumper set
- A paper roll stand and cutter

\*) Multiple material and color options available, please ask for more information

