

Bedside cabinet 2020

Manufactured from high quality materials

The frame of the cabinet table is white epoxy-coated steel. Front panel of the drawers and surface tops are solid birch or beech Iki high pressure laminate. The height of the meal tray is easy to adjust and it can be changed for left-handed and right-handed use without any tools.

The cabinet has four lockable wheels, a safety limiter to avoid the risk of the drawers being pulled out and a slipper rack as standard.



Technical data and features	
Width	49 cm
Depth	62 cm
Height	90 cm
Tray height adjustment range	67–107 cm
Weight	35 kg
Drawer dimensions	28 × 37 × 12 cm
Meal tray	36 × 49 cm
Warranty	2 years
Country of manufacture	Finland

Factory installed options:

- Lock for the top drawer
- Meal tray that can be tilted for use as a reading support
- Chrome-plated side basket and attachment brackets
- Nurse calling device adapter attachment

