

Afia 4060 Gynaecological Examination Chair

Versatile and easy-to-use adjustments



Main features:

1. Motorised adjustment of height, back section and seat section
2. Compact footprint: overall length 140 cm
3. Electronic adjustment to Trendelenburg position

The Afia 4060 is an easy-to-use examination chair developed specifically for gynaecological and urological procedures. With an overall length of 140 cm, it has a small footprint and fits examination rooms where space is limited.

Height, back section and seat section adjust electrically with the hand control, and the 61–98 cm height range gives an ergonomic working position for the clinician. The chair can also be tilted electrically to the Trendelenburg position. A two-part top, adjustable heel supports, a gynaecological bowl, a neck cushion and a paper roll holder are included as standard, and the chair is available in several upholstery materials and colours.

Technical data and features	
Height adjustment	61–98 cm
Patient surface width	59 cm
Patient surface length	122 cm
Overall length	140 cm
Back section length	87 cm
Seat section length	31 cm
Back section adjustment	-7 to +70°
Seat section adjustment	0–25°
Electric Trendelenburg adjustment	0–7°
Total weight	75 kg
Safe working load (SWL)	150 kg
Castors	2 × 65 mm, 2 × 150 mm
Electric system	24V
Other standard features	Hand control; adjustable heel supports; gynaecological bowl, including fixings and adjustment parts; paper roll holder; neck cushion
Upholstery	Future, Nature or Skai Conductive; several colour options
CE marked	Yes (The device is a CE marked medical device according to MDR 2017/745 class 1)
Country of manufacture	Finland

Factory-installed options:

- Arm supports
- Foot controller for height adjustment
- Calf supports

Other accessories:

- Cutter for paper roll holder

