

Capre F Treatment table

Versatility for physiotherapists



Capre F5 table pictured with accessories



Main features:

1. Very low minimum height: Electric height adjustment 45 - 92 cm or hydraulic adjustment 49 - 90 cm
2. Stepless arm supports featuring quick adjustment with gas springs (only F5), anatomic upholstery support the arm of the patient
3. Three or five sections

The modifiability, adjustability and silent operation of the sturdy Lojer Capre F improve functionality and enhance the comfort of use. Its multiple height settings, easy-to-use concave armrests and side supports, and the stepless adjustable extension and flexion angle of the foot section improve the ergonomics and efficiency of the work being carried out.


All Lojer Capre treatment tables have a wide range of upholstery colors and several high-quality upholstery materials to choose from. It is easy to keep the table clean, thanks to its straight surfaces, encased motor and high quality, easy to clean material.

LOJER® *For easy care*

Lojer Capre F	F3	F5	F3H	F5H
Sections	3	5	3	5
Total width	55, 60, 65 and 75 cm (option)	55, 60, 65 and 75 cm (option)	55, 60, 65 and 75 cm (option)	55, 60, 65 and 75 cm (option)
Total length	202 cm	202 cm	202 cm	202 cm
Total weight	86 kg	95 kg	100 kg	109 kg
Height adjustment	45 - 92 cm	45 - 92 cm	49 - 90 cm	49 - 90 cm
Foot section adjustment range	+70°...-12° °	+70°...-12°	+70°...-12°	+70°...-12°
Head support adjustment range	+35°...-90°	+35°...-90°	+35°...-90°	+35°...-90°
Arm support adjustment range	-	34 cm	-	34 cm
Arm support angle adjustment	-	50°	-	50°
Standard features	Stepless head section with gas spring adjustment Head section, oval or U-shape Foot section extension and flexion Electric height adjustment with foot control Foot adjustment (levelling)	Stepless head section with gas spring adjustment Head section, oval or U-shape Foot section extension and flexion Electric height adjustment with foot control Foot adjustment (levelling) Fully adjustable concave arm rests (gas springs)	Stepless head section with gas spring adjustment Head section, oval or U-shape Foot section extension and flexion Hydraulic height adjustment with foot pedal Foot adjustment (levelling) Fully adjustable concave arm rests (gas springs)	Stepless head section with gas spring adjustment Head section, oval or U-shape Foot section extension and flexion Hydraulic height adjustment with foot pedal Foot adjustment (levelling) Fully adjustable concave arm rests (gas springs)
Upholstery	Multiple upholstery colors and several high-quality upholstery materials available	Multiple upholstery colors and several high-quality upholstery materials available	Multiple upholstery colors and several high-quality upholstery materials available	Multiple upholstery colors and several high-quality upholstery materials available
Safe Working Load (SWL)	210 kg	210 kg	210 kg	210 kg
Electrical system	24V	24V	24V	24V
CE-marked according to MDR	Yes	Yes	Yes	Yes
Warranty (years)	2/10	2/10	2/10	2/10
Country of Manufacture	Finland	Finland	Finland	Finland

Accessories:

- Side supports with fast locking and horizontal adjustment in three positions
- Electrical table top height adjustment with 360° foot bar (including safety switch)
- Centrally locking castors (Ø75 mm)
- Double hinged head support
- Drainage with gas spring
- Drainage with electric adjustment
- Additional hand control or foot pedal
- Face hole cover
- Paper roll holder (50 or 60 cm)
- Battery back-up
- Safety switch
- Turning arm rests (can be used as arm and side rests, basic arm rests not available with this feature)
- Extra width upholstery 75 cm
- Conical upholstery shape (only with model F5 and 65 cm width)
- 2-colour upholstery (same as in FX5, only available with F5 and 55, 65, 75 cm widths)
- Ankle cushion



LOJER® For easy care