



Designed for Efficient Patient Care



Product pictured with optional accessories



Main features:

1. Intuitive Basic functions – simple and quick to operate
2. Effortless Maneuverability – smooth transport even in confined spaces
3. Modular and customizable design – flexible configurations and a wide range of accessories to meet diverse clinical requirements

The Emergo patient stretcher is designed to support efficient patient transport and clinical workflows in demanding healthcare environments. Its excellent manoeuvrability allows smooth movement even in confined spaces, while intuitive controls and single-handed side rail operation enable quick and effortless use in everyday care situations.

The open design provides excellent access to the patient from all sides, supporting efficient clinical procedures and ergonomic working positions for healthcare professionals. Easy-to-clean surfaces and a low-maintenance structure help meet high hygiene standards and simplify daily operations.

The modular design, combined with a comprehensive range of accessories, allows the stretcher to be configured for a variety of clinical applications and user requirements. The robust construction ensures durability for intensive use while maintaining ease of handling and patient comfort.

Technical data and features	
Mattress base section	2- or 3-sectional
Width	80 cm
Patient surface width	70 cm
Length	219 cm
Recommended mattress size	70 x 200 cm
Height adjustment	50–87 cm (56–93 cm with X-ray platform)
Safe Working Load (SWL)	250 kg
Trendelenburg	14°
Anti-Trendelenburg	7°
Back section adjustment	0–75°
Leg section adjustment	0–35°
Castors	Central locking, antistatic 150 mm
Country of manufacture	Finland
CE marked	Yes (The device is a CE marked medical device according to MDR 2017/745 class 1)
Warranty	2 years

Factory installed options:

- X-ray platform
- Oxygen bottle holder under the mattress base (2 kg or 5 kg)
- Paper roll holder
- Central locking castors 200 mm XSP
- 5th wheel
- Directional castor (under back or foot section)
- ABS storage box (large)

Other options:

- Push bars: Foldable push bar (pair), foldable wide push bar, wide push bar
- Mattress (thickness 8 or 10 cm)
- IV-pole
- Foldable IV-pole (head-end of the trolley)
- Restraining straps
- Infusion rod holder (attached on the side of foot section)
- Accessory basket
- Wire basket (under back/foot section)

